






Rajdhani Surgicals


VETERINARY INSTRUMENTS & EQUIPMENTS

CATALOGUE





Rajdhani Surgicals



rajdhansurgical30@gmail.com
rajdhani_surgical@yahoo.com

Anesthesia Machine for animal



1001

Features & Benefits:

- Simple. Effective. Easy.
- 0 – 4 LPM flowmeter
- Free Darvall ZDS Mask non-rebreathing circuit (#9130)
- Oxygen supply hose
- Instruction manual
- Responsive technical support
- Table-top / Wall Mount model
- Increase mobility by adding a compact 5-Star Wheelbase (#1569)
- DVM Iso Vaporizer (#9176) sold separately

Color Doppler Ultrasound



Awarded by Govt. of India



Price on demand

Ethioli Tiny 16a Classic Trolley with LED Monitor

1002

- Colour Doppler, Pulse Wave, B, B/M Imaging modes
- Large 15" LED flat panel monitor with tilt and swivel for Comfortable Viewing
- 30cm wide ergonomic and intuitive keyboard
- Report entire scan on CD/DVD as a video
- Assistant free reporting - hand out structured reports with images right at the end of the exam
- Multi frequency broad band probes; Tissue Harmonic
- Imaging; Speckle reduction imaging
- Three probe port connectivity

Ultrasound Portable Model



1003

Electric Animal Operation Table




1004

SPECIFICATIONS

- Electric Animal Operating Table • Electric up and down function.
- It is suitable for small animal's operation, such as sheep, dog etc.
- It is reasonable structured, tight, high performance, easy operation and has very strong practicability.
- Electric animal operating table is composed of table-board, pin, and table pedestal, equipped with surgical instruments bracket, tray and infusion support.
- The completed operating table is made of stainless steel, table pedestal adopts rolling design, can move easily and floor locked.
- Fully made of Premium Grade 304 Stainless Steel, it is heat-resistant, corrosion protective and rustproof.
- The surface of the table-top is V shaped and is smooth and easy to clean and disinfect, thus, it can avoid cross infection • Equipped with fluid collecting through system.

Surgical Operating Table Hydraulic Side end control



1005

Small Animal
Hydraulic Table
Side & Control

SALIENT FEATURES

Small animal Veterinary surgical table, The table can be rotated 360 degrees Can be tilted 45 degrees, front and back stainless
BS0046 Small animal Veterinary surgical table
1. The table board is U-type, and uses hydraulic system to control rise/down, equipped with drainage port. 2. The table can be rotated 360 degrees. 3. The working table can tilted about 15 degrees, left and right, mechanic manual operation. 4. The table can be tilted 45 degrees, front and back, respectively, manual operation. The completed table is made of stainless steel. 5. The pedestal is equipped with a pulley.
Table board dimension (L X W) 1200 X 500 mm Table height : 820-930mm

Large Animal Manual Table with 2 wheels & 2 Handle



1006

Specification

Top SS
4 Hooks at Back Side
Dimension (L X W) 8ft X 6ft
Table height : 1.75ft

Electric Animal Operation Table



1007

Small Animal Examination Table (top Ss Base Powder Coated With Rubber Matt)



1008

SPECIFICATION :

Top : S.S. Grade : 304
L 4' x W : 2' x H : 3'

Sheet : 20 Gauge
Base Frame : 18 Gauge

Powder Coating
Wash Basin Provision

I.V. Attachment
Inner Side Hooks

Visual Ai Gun

**For performing
Artificial Insemination
without Rectal
manipulation with
LED Light / Camera**



MADE IN GERMANY / INDIAN

1009

LED Lights



Specifications

Intensity	: 90,000 + 90,000 Lux \pm 10%
Size of Light Field	: 12-30 cm
Colour Temperature	: 4000-5000 K
Colour Reduction Index	: 93 RA
LED Life	: >50000 Hrs.
Number of LED	: 84+63
Diameter of Light	: 500mm
Brightness Control	: Capacitive Touch Panel
Power Supply	: 220v/50Hz AC
Focusing	: Adjustable
Battery Backup	: Optional

1010



Mobile Surgical Light With Battery Backup

Specifications

Intensity	: 90,000 Lux \pm 10%
Size of Light Field	: 12-30 cm
Colour Temperature	: 4000-5000 K
Colour Reduction Index	: 93 RA
LED Life	: >50000 Hrs.
Number of LED	: 63
Diameter of Light	: 500mm
Brightness Control	: Capacitive Touch Panel
Power Supply	: 220v/50Hz AC
Focusing	: Adjustable
Battery Backup	: Optional

1011



Specifications

Intensity	: 90,000 Lux \pm 10%
Size of Light Field	: 12-30 cm
Colour Temperature	: 4000-5000 K
Colour Reduction Index	: 93 RA
LED Life	: >50000 Hrs.
Number of LED	: 63
Diameter of Light	: 500mm
Brightness Control	: Capacitive Touch Panel
Power Supply	: 220v/50Hz AC
Focusing	: Adjustable
Battery Backup	: Optional

1012



Specifications

Intensity	: 90,000 Lux \pm 10%
Size of Light Field	: 12-30 cm
Colour Temperature	: 4000-5000 K
Colour Reduction Index	: 93 RA
LED Life	: >50000 Hrs.
Number of LED	: 63
Diameter of Light	: 500mm
Brightness Control	: Capacitive Touch Panel
Power Supply	: 220v/50Hz AC
Focusing	: Adjustable
Battery Backup	: Optional

1013

Easily detachable connector based OT Light / Optional



Specifications

Dome Size	: 500 mm
Intensity Lux	: 140000 \pm 10%
Spot Diameter	: 125-200mm
Colour Temp.	: 3500k-4800k
Focusing	: Manual
Power Supply	: 220v /50Hz AC
Colour Rendering Index	: 95RA
Power Consumption	: 40Watt
Led Average Light	: 50000Hrs
Intensity Control	: 30~100

1014



Specifications

Dome Size	: 500 mm
Intensity Lux	: 140000 \pm 10%
Spot Diameter	: 125-200mm
Colour Temp.	: 3500k-4800k
Focusing	: Manual
Power Supply	: 220v /50Hz AC
Colour Rendering Index	: 95RA
Power Consumption	: 40Watt
Led Average Light	: 50000Hrs
Intensity Control	: 30~100

1015



Specifications

Dome Size	: 500 mm
Intensity Lux	: 140000 + 140000 \pm 10%
Spot Diameter	: 125-150mm
Colour Temp.	: 3500k-4800k
Focusing	: Manual
Power Supply	: 220v /50Hz AC
Colour Rendering Index	: 95RA
Power Consumption	: 40Watt
Led Average Light	: 50000Hrs
Intensity Control	: 30~100

1016



Specifications

Dome Size	: 500 mm
Intensity Lux	: 140000 \pm 10%
Spot Diameter	: 125-150mm
Colour Temp.	: 3500k-4800k
Focusing	: Manual
Power Supply	: 220v /50Hz AC
Colour Rendering Index	: 95RA
Power Consumption	: 40Watt
Led Average Light	: 50000Hrs
Intensity Control	: 30~100

1017



Specifications

1018

Intensity	: 50,000 Lux \pm 10%
Size of Light Field	: 12-30 cm
Colour Temperature	: 4800 K
Colour Reduction Index	: 93 RA
LED Life	: >50000 Hrs.
Number of LED	: 21
Diameter of Light	: 150mm
Brightness Control	: Touch
Power Supply	: 220v/50Hz AC
Focusing	: Adjustable



Specifications

1019

Intensity	: 45,000 Lux \pm 10%
Size of Light Field	: 12-30 cm
Colour Temperature	: 4800 K
Colour Reduction Index	: 93 RA
LED Life	: >50000 Hrs.
Number of LED	: 14
Diameter of Light	: 150mm
Brightness Control	: Touch
Power Supply	: 220v/50Hz AC



Specifications

1020

Intensity	: 40,000 Lux \pm 10%
Size of Light Field	: 12-30 cm
Colour Temperature	: 4800 K
Colour Reduction Index	: 93 RA
LED Life	: >50000 Hrs.
Number of LED	: 14
Diameter of Light	: 150mm
Brightness Control	: Capacitive Touch Panel
Power Supply	: 220v/50Hz AC
Focusing	: Adjustable

Led Light Source



1021



Halogen Light Source



1022

Ambu's Respirator for veterinary Set

Powerful 120W endoscopy LED light source with the latest LED technology . Equivalent lighting to 350W Xenon light source. Perfect for use with all types of endoscopes - rigid and flexible. Remote control is provided for the ease of use, able the doctor to adjust the setting without leaving the patient chair.

LED Head light Rechargeable wireless



1023



1024

Small Pet / Large Feline
32490B16
2.5" H x 3.75" W
Diaphragm Opening 1.125"



Medium Pet / Small Canine
32490B15
3" H x 3.5" W
Diaphragm Opening 1.125"



Large Pet / Large Canine
32490B10
4.5" H x 5.5" W
Diaphragm Opening 2.25"



Large Animal Mouth Gag

Material
Copper
Gun Metal
Nickel
Silver 100gms
Weight
4.100kg
Approx

1025



Horse Mouth Gag

MADE IN GERMANY
1026



Weingarth Rumenotomy Set



1027

Stainless steel.

Ob Chain Stainless 60"



1028

MADE IN GERMANY / INDIAN

Ob Chain Handle



1029

Stainless steel.
4" long x 1/2" diameter.
MADE IN GERMANY / INDIAN

Krey Hook



MADE IN GERMANY / INDIAN 1030

Mouth Gag For Cat & Dog



MADE IN GERMANY / INDIAN


1031

Harm's OB Eye Hook Blunt / Sharp



1032

OB Wire Handle



Sutter Style OB Handle

1033

OB WIRE GUIDE

Hauptner OB Wire Guide



1034

OB Wire Drop Guide



1035

OB Wire Guide



1036

OB WIRE & HANDLE

Hauptner OB Wire



1037

Buhner Suture Needle



1038

Veterinary Stomach Pump



Indian / Imported
1039

ET Tubes



1040

Ai Gun Universal



With Lock &
Without Lock


1041

Ai Sheath



1042

'Canula Needle



1043

'Centasis needle'

(Reducer)
for farm animals
pericardiocentesis,
abdomincentesis,
Renal and liver biopsy.

Re-usable
1044



Urinary Inconsistence Kit



1045

URINARY INCONSISTENCE KIT CONSISTING OF :
1. SILICONISED CATHETER WITH INTRODUCER
2. URINE COLLECTING BAG
Size Available : 8 FR, 10 FR, 12 FR, 14 FR, 16 FR

Suture Needles



All Sizes Available Round Body & Cutting
Packet of 6 Needles

1046

Infrared Lamp



-Used for Heat Therapy
-Treats Muscular and Joint Pain

1047


Teat Surgery Kit



Stainless steel.

1048

X-ray view box



1049

Nasogastro Tube for animal



Indian / Imported Available **1050**

Stomach Tube for animal



Indian / Imported Available **1051**

Gyn Stick



MADE IN GERMANY

Product Description

A new OB multi-use instrument designed to address a number of bovine dystocia issues. The Gyn-Stick can help as

– Uterine torsion correction – Kuhn's crutch for repelling a malpositioned calf. Simply turn the Gyn-Stick around, insert rectally in breech calf. Instructional video at www.JorVet.com.4

1052

Bovine Cervical Forceps



MADE IN GERMANY / INDIAN

1053

Castration



9" & 14"

1054

Animal Body Metal Detector



Alarm Mode : Sound
Features : Rechargeable Batteries
Range : 2 meters

Made From MS Body Electric
Use in Animal Body Metal Testing
Available Size : Regular
with 2 Chargers

1055

Trephiner With Handle



1056

Postmortem Of Large Animal



1057

Laryngoscope For Small Animals



1058

Laryngoscope For Large Animals



1059

Freemartin Probe



1060

Thygesen's Fetatome



Product Description

Two-tube model with adjustable handle for OB chain. This is the breakdown model, complete with insertion coil and cleaning brush. 33L.

MADE IN GERMANY / INDIAN

1061

Ultrasonic Doppler Animal Blood Pressure System

The Vet-Dop2 is used to screen for hypertension, to check blood pressure in surgery, to monitor blood flow at the extremities during surgery and to check for intact blood vessels after trauma and before amputation. The non-invasive Doppler technique is the method of choice for measuring blood pressure in animals weighing less than 15 pounds, especially cats. Other methods, including monitors that use the oscillometric technique with small animal validated software may not calculate pressures as quickly especially on awake, smaller patients. The Vet-Dop2 also works well on larger breeds and can detect blood flow at the base of the tail in horses. The Vet-Dop2 is the Doppler of choice and preferred by experts in the field of veterinary blood pressure measurement.



1062



HEINRICH



High Quality

Reliable

*For a Better Life
For our Loyal Friends*

Air Drill



1063

Battery Drill



1064

Tplo Saw



1065



Our Products are Lightweight and Portable. They provide the best solutions in different areas of veterinary medicine.

HEINRICH

MEDIZINTECHNIK Heilbronner strasse 150,
Stuttgart 70191, Germany Website : www.heinrichmedical.de

Veterinary Bone Drill



Battery Operated


Made in Korea




1066

No.	Specification
M-11	<ol style="list-style-type: none"> 1. Maximum speed: 1200rpm, stepless speed regulation, rotating in positive and negative direction for selection, stable operation 2. Gripping drilling bit diameter: 0.4-4mm 3. Battery voltage: 7.2V; free from disinfection; 2 batteries. 4. Reliable performance; lightweight; suitable for hand, foot, frontal plane surgical operation; 5. Strong host power, conjoined closed cover 6. Suitable for high temperature sterilization

Manual Bone Drill



1067



1068

Veterinary Pulse Oximeter with Thermistor



1069

Details

The Model 340V is unique with excellent low perfusion monitoring, sturdy washable clip-on sensor, rectal probe with thermistor for oxygen saturation plus body temperature (in degrees celsius or Fahrenheit), and automatic temperature display every two minutes or when solicited.

Dimensions

7.5 (L) x 3.5 (W) x 1.4 (H) in.
19.1 (L) x 8.9 (W) x 3.5 (H) cm.

Technical Specifications

Display: Two 3-digit, 7-segment LED
Weight: 18.2 oz. (516 grams)

Resolution:

SpO2: 1%
Pulse: 1 BPM
Temperature: 1°C / 0.2°F

Range:

SpO2: 0-100%
Pulse: 32-250 BPM
Temperature: 24.5°C - 42.5°C, 76.1°F - 108.5°F

Accuracy:

SpO2:
100-70% ± 2%
69-60% ± 3%
Less than 60% unspecified
Pulse: 32-250 BPM ± 2 BPM

Temperature:

24.5°C - 42.5°C
76.1°F - 108.5°F

Power Source:

NiCad rechargeable battery - lasts 12 hours in continuous operation.
Power adapter - 100 - 250 V AC 50/60 Hz, 0.5 A
Warranty Information
Five-year limited warranty.

IP300 Veterinary Infusion Pump



1070

Specification

Model Number : IP300
Instrument Classification : CLASS I
Display : Large LCD Display
Packing Size : 300 x 260 x 160 mm
G.W. : 2.5 Kg
Work Modes : : Five Work Modes

ECG300G-Vet Animal Veterinary ECG Machine



1071

Specification

Type : Portable ECG / EKG Animal / Veterinary Use
Syringe : Portable ECG / EKG Animal / Veterinary Use
Blood Sampling Needle : Portable ECG / EKG Animal / Veterinary Use
Breathing Pattern : Portable ECG / EKG Animal / Veterinary Use
Product name : Portable ECG / EKG Animal / Veterinary Use
Channel : : 3

ECG300G-Vet Animal Veterinary ECG Machine



1072



1. Introduction

ECG100G-VET single channel electrocardiograph can print ECG waveform with thermal printing system, prompt for "Lead off" and "Lack of paper", display multi-language interface, record ECG waveform with manual or auto-mode, which is very convenient for operation.

Specification

Type : Portable ECG / EKG Animal / Veterinary Use
Syringe : Portable ECG / EKG Animal / Veterinary Use
Blood Sampling Needle : Portable ECG / EKG Animal / Veterinary Use
Breathing Pattern : Portable ECG / EKG Animal / Veterinary Use
Product name : Portable ECG / EKG Animal / Veterinary Use
Channel : : 3

Ophthalmoscope



Beta - ZUO
HEINEZ
1073

Otoscope



Mini - 2000
HEINEZ
1074



Mini - 2000
HEINEZ
1075

Made in Germany

Otoscope



INDIAN
1076

Otoscope & Ophthalmoscope ENT Set



Small Animal ENT Set IMPORTED

1077

CMS8000VET Medical Equipment Animal/Veterinary Monitor

Introduction

CMS8000VET Medical Equipment Animal/Veterinary Monitor Multi-parameter Monitor has abundant monitoring functions and can be used for monitoring of cat, dog and other animals. User may select different parameter configurations according to different requirements. Veterinary Monitor adopts 100-240V~, 50/60Hz for power supply, 12.1" TFT color LCD for displaying real-time data and waveform, and up to 8-channel waveform and all monitoring parameters can be displayed simultaneously. In addition, Veterinary Multi-parameter Monitor can be connected to the central monitoring system via wired/wireless network so as to form a network monitoring system. Animal Monitor device can monitor such parameters as ECG, RESP, NIBP, SpO2 and dual-channel TEMP, etc. It integrates parameter measurement module, display and recorder in one device to form a compact and portable equipment. At the same time, its built-in replaceable battery provides convenience for moving.

Functions

Standard parameters: ECG, RESP, SpO2, PR, NIBP, dual-channel TEMP

- 1) ECG Heart rate (HR)
ECG waveform
Arrhythmia and ST-segment analysis
- 2) RESP Respiration rate (RR)
Respiration waveform
- 3) SpO2 Oxygen saturation (SpO2)
Plethysmogram (PLETH) waveform
Pulse rate (PR)
Bar graph
- 4) NIBP Systolic pressure (SYS), Diastolic pressure (DIA), Mean pressure (MEAN)
- 5) TEMP T1, T2, TD

It has abundant functions, such as audible and visual alarm, trend data storage and output, NIBP measurement, alarm event marking and drug concentration calculation, etc

Features

- 1) 12.1" TFT color LCD, multi-language interface (Simplified Chinese, Traditional Chinese, English, French, German, Turkish, Spanish, Portuguese, Italian, Dutch, Romanian, Russian, Kazakhstan, Polish, Czech).
- 2) Fanless design, quiet, energy-saving and clean, which reduces the possibility of cross-infection.
- 3) With standard interface, OxyCRG, trend chart, big character interface and view bed, convenient to observe.
- 4) Operation mode: keys and buttons.
- 5) Display up to 8-channel waveform on one screen.
- 6) Display 7-lead ECG waveform on one screen, and with ECG step function.
- 7) Adopt digital SpO2 technology, anti-motion and anti-ambient light interference, and measurement can be performed under the circumstance of weak filling.
- 8) NIBP measurement mode: Manual/AUTO/STAT, storage for 4800-group NIBP data.
- 9) Review for 71 alarm events of all parameters and 60 arrhythmia alarm events.
- 10) Drug concentration calculation and titration table functions.
- 11) Print trend table data with one-key.
- 12) Connect to Central Monitoring System by 3G, WiFi or wired mode.
- 13) AC/DC, built-in rechargeable lithium battery achieve uninterrupted monitoring.
- 14) Anti-high frequency surgical unit, defibrillation-proof (special leads are necessary).



1078

Pet Blood Pressure Monitor 2.8" color LCD screen -CONTEC08A-VET

CONTEC



1.Introduction

Pet Blood Pressure Monitor 2.8" color LCD screen -CONTEC08A-VET can measure blood pressure for animals **with 5 cuffs**.

This Clinic veterinary instrument features in full-automatic measurement, high-definition color LCD display and strong visibility.


- with integrated parameter measure, display and record output in one device,
- CONTEC08A-VET pet blood pressure monitor Adopts data review interface such as "data list", "trend chart", "big font" with easy and simple operation.
- be applied to routine test of family pets, animal hospital, clinic and farm, etc.
- Spo2 probe is optional

2.Features

1. Three cuff mode: large, medium, small. User can choose corresponding cuff and cuff mode according to different body type of animals.
2. High-definition color LCD display, strong visibility.
3. Data storage function, up to 100 groups of data can be stored, by which measured data can be reviewed.
4. Physiological alarm function, alarm limits can be set the alarm limits. When blood pressure is higher than the high limit or lower than the low limit, the alarm will occur. Alarm switch can be set.
5. Supply two kinds of units: Kpa/mmHg, which can be set as user's need. Audial and visual prompt when power is low.
6. With short-time power storage function, when replacing the battery, the clock time function can't be affected.
7. One-touch operation, automatic BP measurement, adopts oscillometric multistep deflation measurement method, good repeatability. When measurement error happens, device will prompt error message.
8. Store measured results with date and time.
9. Function of automatic power-off.

1079

Microcontrolled Electrosurgical Diathermy-M 250



1080

Micro Controlled Surgical Diathermy Unit



DELCAIT-MP-4

1083

Micro controlled Electro Surgical Diathermy with Bipolar TURP "elian"



1084

Surgical Diathermy Model :

Delcatt - 250B
DELTA MEDICAL CORPORATION
 Manufacturing Electro-Surgical and Physiotherapy Equipments



1081

CO2 Insufflicator



1085

Vessel Sealer With Micro Controlled Surgical Diathermy



1082

Cervical Cautery



1086

4Ltr Capacitor Ulv Fogger



Features

Diameter : 23.5 cm
Capacity : 4 Litre
Motor Speed : 20000 rpm
Material Type : Body & Tank: 304 Grade Stainless Steel
Dimensions : 33 x 41 x 23.5 cm

1087

Shoe Despensor



Easy to carry
Prevents contamination
with 200 pcs. Shoe Covers

1088

Stethoscope



Indian / Imported Available

1089

Teaching Stethoscope



Indian / Imported Available

1090

Fluid Suction Machines



HI VAC SUCTION MACHINE

1091



SLOW SUCTION MACHINE

1092



SUCTION MACHINE

1093

TRAY LIFTING STERILIZER 20"



1094

DRESSING DRUM ISI MARK 11X9



1095

BUCKET



1096

CHETTLER JAR



1097

**AUTOClave SINGLE DRUM
STEEL FINISH ALUMINUM**



1098

**AUTOClave DOUBLE DRUM
STEEL FINISH ALUMINUM**



1099

SS TRAY 10 X 12



1100

KIDNEY TRAY



1103

**AUTOClave DOUBLE DRUM
STEEL FINISH (ALL SS)**



1101

**ELECTRIC STERILIZER
SIZES AVAIL : 10", 12", 14", 20"**



1102

SS TRAY 8 X 10



1104

Mayos Trolley (All S.S.)



1105

Instruments Trolley (All S.S.)



Size : 18X24

1106

Wash Basin Stand (all S.s.)



1107

I.V. STAND



1108

Specification

L x B x H	50 x 6 x 6 inches
Weight	5Kg

Salient Features

- Stainless steel pole
- Heavy duty ABS portable base
- 3" smooth glide casters
- 4 U-hooks

STOOL 4 LEGS



STOOL 3 LEGS



STOOL WITH CUSHION



4 FOLD SCREEN



MEDICINE TROLLEY



Leg Band



1114

Wing Band




1115

Aluminum Ear Tag



1116

Egg Tray



1117

Automatic Bell Drinker



1118

Rebreathing Bag



15 / 30 Ltr. 1119
IMPORTED MADE IN USA

Automatic Drinker



1120

Mini Automatic Drinker



1121

Poultry Grower Feeder



1122

Poultry Chick Feeder



1123

Poultry Cage Drinker



1124

Poultry Chicks Drinker



1125

Grower Drinker



1126

Poultry Plastic Fogger



1127

Sprinkler System



1128

Poultry Automatic Vaccinator




1129

Flame Gun



1130

Poultry De-Beaker Machine



1131

Pulley



1132

Automatic Drencher



1133

Imported
Available in all sizes

Poultry Defeathering Machines 5 to 7 Birds



1134

Grain Moisture Meter



1135

Infrared Thermometer SRG 900



1136

TDS/Conductivity/Temp. Combo Meter Thermometer



1137

Electronic Distance Measuring Tool



1138

Dog & Cat Refractometer



1139

Electric Brooder



1140



Formalin Chamber

(SIZE 14", 20", 23", 26")
6mm Sheet Heavy Duty

1141

Needle Syringe Destroyer (SUPER DELUXE)



1142

Digital Conductivity Meter



Specifications	
Type	Digital Conductivity Meter
Certification	CE Certified
Weight	0.8 Kg approx.
Color	White
Voltage	3-5VDC, 6-9VDC
Usage	Industrial Usage
Display Type	Digital
Automatic Grade	Automatic
Power Source	Battery, Electric
Features	Accuracy, Durable
Warranty	1Year
Display	3 1/2 Digit LED
Power	12V DC using 230V + 10% AC Adapter
Dimensions	270x202x80 mm (bxbxh)
Cell Constant	Adjustable on Digital Display

1143

Small Animal Instrumental (Kit of 45 Pcs)

Our Small Animal Kit Consists of:

- 1 x Stainless Steel Tray with Lid
- 1 x Scalpel Handle No 3
- 2 Dissecting Forceps 6"
- 2 Dissecting Forceps 4"
- 4 x Allis Tissue Forceps 3/4 Teeth 15cm 6"
- 2 x Needle Holder 14cm 5.5"
- 1 x Mayo Scissors 14cm 5.5" Str
- 4 x Towel forceps 4"
- 4 x Towel Clips cross sections
- 4 x Halstead Mosquito Forceps 13cm 5" Cvd
- 2 x Rochester Pean Artery Forceps 16cm 6" Cvd
- 2 x Skin Hook Single
- 2 x Single Hook Double
- 1 x Sponge Forceps 8"
- 1 x Aneurysm needle
- 2 x Babcock forceps 6"
- 4 x Kocher forceps 6"
- 2 x Trochar Canula
- 4 x IS Clamp Ped.

1144

Large Animal Instrumental Set

Set of 105 Instruments

- 2 pcs Surgical scissors (curved /round), 8"
- 2 pcs surgical scissors (straight /round), 8"
- 4 pcs surgical scissors (straight /curved), 12"
- 6 pcs hemostatic forceps (straight, curved),
- 4 pcs dissecting forceps 6"
- 4 pcs dissecting forceps 8"
- 2 pcs dissecting forceps 12"
- 4 pcs Skin hook single
- 4 pcs Skin hook double
- 2 Pc Needle holder 6"
- 2 Pcs needle holder 8"
- 2 Pcs needle holder 10"
- 4 pcs Allis forceps 6"
- 4 pcs Allis forceps 8"
- 4 pcs forceps 8"
- 2 pcs Babcock forceps 8"
- 2 pcs Babcock forceps 10"
- 4 pcs Kocher forceps 6"
- 4 pcs Kocher forceps 8"
- 4 pcs I S Clamp 10"
- 1 pcs Ribbon retractor maliable,
- 4 pcs Artery Forceps 6"
- 4 pcs Artery Forceps 8"
- 4 pcs Artery Forceps 10"
- 3 pcs Devers retractor
- 2 Doyan retractor
- 2 pcs Morries retractor
- 4 pcs Right Angle retractor
- 2 pcs London Hospital Angle retractor long blade
- 1 maliable retractor
- 6 towel forceps
- 2 Scalpel handle
- 4 Bone trephine
- 4 pcs Sponge holding forceps
- 1 SS Tray

1145

Large Animal Orthopaedic Kit (Set of 90 pcs)

- 1 pc Electric bone cutter wit blades
- 1 pc Electric bone drill
- 2 pcs Arthroplasty Rasp
- 1 pc Universal Rasp 210mm Long
- 2 pcs Hohman Retractor 18mm Short Narrow Tip
- 2 pcs Hohman Retractor 12mm Short Narrow Tip
- 1 pc Standard Jacobs Chuck & Handle 1/4" (6mm) Capacity
- 1 pc Bone Cutter S/Steel Cuts To 4.0mm (370mm Long) TC Inserts
- 1 pc Bone Cutter S/Steel Simple Action 2.0mm Max (255mm Long)
- 1 pc Wire Passer Double Ended (145mm)
- 1 pc Wire Twister & Shear Cutter (165mm)
- 1 pc Periosteal Elevator CVD AO Type 180mm 6mm End PTFE Handle
- 1 pc Small Curved Rongeurs 4mm Bite
- 1 pc Small Angled Cutter (15mm Blade)
- 1 pc Large Angled Cutters 19Mmm Blade 23cm Long
- 1 pc Bone nibbler straight & curved
- 1 pc Screw box
- 1 pc K nails set of 10 pcs
- 1 pc Mallet 300g
- 1 pc Mallet 140g
- 1 pc 10mm Osteotome Light (140mm)
- 1 pc Kirschner & Arthrodesis Wire Bender 140mm
- 1 pc Stubby Punch 50mm
- 1 pc Modular Osteotome With Seven Blades - 4,6,8,10,12,15 & 20mm
- 1 pc 1.1mm Arthro Wire x packet 10
- 20 pcs K/A Wire
- 1 pc 15mm Osteotome Swedish Pattern (205mm)
- 1 pc Volkman Scoop Double Ended 4mm & 6mm (200mm Long)
- 20 pcs Drill bits of different sized
- 30 pcs Steinman Pins Trochar Tips
- 10 pcs Reel Of Orthopaedic Wire Diameter

1146

Large / small Animal Ophthalmology Kit (Set of 38 pcs)

- 4 pcs Halstead Mosquito Forceps 13cm 5" Str
- 4 pcs Halstead Mosquito Forceps 13cm 5" Cvd
- 1 pc Lims Forceps
- 1 pc Eyelid Speculum (bigger size spec for larger animals)
- 2 pc Wire Speculum large / small
- 1 pc Castroviejo Needle Holder 14.5cm 5.5"
- 2 pcs Corneal scissors, 11cm Right / Left
- 1 pc Corneal scissors, 11cm Straight 5"
- 1 pc Corneal scissors, 11cm Curved 5"
- 2 pcs Suture tying forceps straight / curved
- 1 pc Mech pherson forceps
- 1 pcs Iris repositer
- 1 pc Epyliation forceps
- 1 pc Iris forceps curved
- 2 pcs vannus scissors straight / curved
- 1 pc fixation forceps
- 1 pc pigtel prob
- 1 pc Puntum dialator
- 1 pc Lacrimal prob set
- 1 pc Lense hook
- 1 pc Eye vectus
- 2 pcs simco canula
- 2 pcs air injection canula
- 1 pc I ball coutary
- 1 pc Wett filled coutary
- 1 pc silicon tray with matt

1147

Small Animal Orthopaedic Kit

No.	Name	Qty.	No.	Name	Qty.
1.	AO A.C. drill bit 1. 1*70	2	31.	small reduction forceps with points	1
2.	AO A.C. drill bit 1. 5*70	2	32.	middle reduction forceps with points	1
3.	AO A.C. drill bit 1. 8*120	2	33.	small reduction forceps curved	1
4.	AO A.C. drill bit 2. 0*100	2	34.	needle reduction forceps curved	1
5.	AO A.C. drill bit 2. 4*120	1	35.	self centering bone holding forceps	2
6.	AO A.C. drill bit 2. 5*110	1	36.	universal palte bender	1
7.	AO A.C. drill bit 3. 5*110	1	37.	4 holes drill guide 2.0	1
8.	AO A.C. cortial screw tapl. 1.5	1	38.	1.5/2.0 double drill guide	1
9.	AO A.C. cortial screw tapl. 2.0	1	39.	1.8/2.4 compression double drill guide	1
10.	AO A.C. cortial screw tapl. 2.4	1	40.	2.0/2.7 compression double drill guide	1
11.	AO A.C. cortial screw tapl. 2.7	1	41.	2.5/3.5 compression double drill guide	1
12.	AO A.C. cortial screw tapl. 3.5	1	42.	screw holding forceps	1
13.	AO A.C. cortial screw tapl. 4.0	1	43.	osteotome	1
14.	AO Q.C. four-angle 1.2 screw driver	1	44.	bone currette	1
15.	AO Q.C. four-angle 1.5 screw driver	1	45.	probe	1
16.	AO Q.C. Hex 2.5 screw driver	1	46.	depth gauge with kirschner wire	1
17.	AO Q.C. countersink drill 5mm	1	47.	ruler	1
18.	AO Q.C. countersink drill 6mm	1	48.	double joint rongeur	1
19.	AO Q.C. screw extraction 4.0	1	49.	thred guide	1
20.	AO Q.C. round head straight handle	1	50.	nipper pliers	1
21.	AO Q.C. T-shape adapter	1	51.	wire & plate cutter	1
22.	2.5mm drill guide	1	52.	16cm skin retractor	2
23.	hammer 180g vood handle	1	53.	17cm L-type skin retractor	2
24.	periosteum elevator 6.5	1	54.	mounth gag	1
25.	periosteum elevator 3.0	1	55.	kirschner wire bender	1
26.	bone elevator 15mm	1	56.	bone file (flat)	1
27.	bone elevator 8.5mm	1	57.	double head retractor	1
28.	deepth gauge	1	58.	wire pliers	2
29.	plate holding forceps	1	59.	template-short (blue)	1
30.	Mini reduction forceps with points	1	60.	reconstructions template (black)	1
			61.	sterilizing case	1

1148

Mini Multifunction Surgical Power Tool

CODE : 1149

Mini Sagittal Saw Attachment
EVM1-S010S

Mini AO Coupling Drill Attachment
EVM1-A010S

Mini Bone Drill Attachment
EVM-020S

K-wire Attachment
EVM-K01S



EV-MS502

Specification

Electric Handpiece:

Weight (with battery): ≤ 500Grams
Speed: 0-16000rpm
Output power: 80W
Batter (Ni-Mh):
Working Voltage: 14.4V
Charging time: 2-3 hrs
Run time with battery: 40-50 min

Specification

Mini sagittal saw attachment

Speed/Oscillations: 0-1600 osc/min
Clamping range:
Working Voltage: 14.4V
Sterilization: up to 135°C
Noise: ≤50dB

Specification

K-wire attachment

Speed : 0-1100 rpm
Working Voltage: 0.8-1.8mm
Sterilization: up to 135°C
Noise: ≤30dB

Features:

1. Applied to hands & feet surgery.
2. Light handpiece, easy carrying multifunction and quick-coupling.

Specification

Mini bone drill attachment:

Speed: 0-1100rpm
Working Voltage: 14.4V
Sterilization: up to 135°C
Noise: ≤30dB

Specification

Mini AO coupling drill attachment

Speed : 0-1100 rpm
Working Voltage: 14.4V
Sterilization: up to 135°C
Noise: ≤30dB



Battery



Charger



Battery Aseptic Trabsery Ring




Key For Drill Chuck



Aluminum Case




Sterilization Case



Sternum saw blade



Reciprocating saw single face tooth



Reciprocating saw double-sided tooth



Rajdhani Surgicals

Manufacturer, Supplier, Importer & Exporter of :

- DERMATOLOGY INSTRUMENTS • PLASTIC SURGERY • PAEDIATRIC SURGERY
- GYNAECOLOGY • E.N.T. • OPHTHALMOLOGY • NEURO SURGERY
- CARDIO THORACIS • LAPROSCOPY • VETRINARY INSTRUMENTS
- NURSERY EQUIPMENTS • PHYSIOTHERPY EQUIPMENTS & SURGICAL COUTARY
- HOSPITAL FURNITURE & O.T. HOSPITAL EQUIPMENTS ETC.

Head Office : G 222, First Floor, Sec-5, Samosa Chowk, Bawana Indl. Area, Delhi

Patna Office : Dahibhatta Kothi, East of Old Agam Kuan, Thana B.H. Colony,
Kankarbagh, Patna, Bihar-800026

E-mail Id : rajdhanisurgical30@gmail.com

T : 011-45056248

M : +91-9312234980

M : +91-9911972372