

DELHI - GSTIN No : 07ADHPC7238 E1ZS  
PATNA - GSTIN No : 10ADHPC7238E1Z5



# Rajdhani SURGICALS

An ISO 9001:2008 Certified Company  
An ISO 13485:2012 Certified Company

**IMPORTER & EXPORTERS OF SURGICALS INSTRUMENTS & EQUIPMENTS**

## URETHROPLASTY COMPLETE SET



Please feel free to Contact us :

Head Office : G 222, First Floor, Sec-5, Samosa Chowk, Bawana Indl. Area, Delhi

Patna Office : Dahibhatta Kothi, East of Old Agam Kuan, Thana B.H. Colony,  
Kankarbagh, Patna, Bihar-800026

E-mail Id : rajdhanisurgical30@gmail.com

T : 011-45056248


M : +91-9312234980

M : +91-9911972372

## Urethroplasty Complete Set

		Qty.
1	Needle holder 3-0 6 inch T.C. Gold	1
2	Ryder Needle holder 4-0 8 inch	1
3	Needle holder Curved 8 inch	1
4	Mastoid Retractor	2
5	Metzenbaum Scissor Curve 8 inch	2
6	Fine Tooth Forceps	2
7	Fine Plain Forceps	2
8	Mosquito Forceps (5") with 1 - 6 marked for easy identification for suturing	10
9	Kristeller Speculum	1
10	Nasal Speculum Fiber Optic	1
11	Turner Warwick Ring Retractor Set with Silicon Skin Hooks	1
12	C Retractor Set of 3	1
13	Boyles Davis Mouth Gag (Set)	1
14	Urithal Dilator Set	1
15	Hey Groove Urethral Dilators Set	1
16	Hammer	1
17	Cheesle with wooden handle	1
18	Silicone Skin Hooks	6
19	Perforated Box	1

## Turner Warwick Ring Retractor with 8 blades



**C Retractor  
Set of 3**



**Kristeller  
Speculum**



**Nasal Speculum  
Fiber Optic**



**Mastoid Retractor**



### Urethral Dilator Set

**Specifications :**

Dilators - 8, 10, 12 and 14 Fr, with sideports, 37 cm long, polyethylene

Dilators - 16, 18, 20, 22, and 24 Fr, without sideports, 37 cm long, polyethylene

Straight safety wire guide - .038 inch (.97 mm) diameter, stainless steel, 80 cm long wire guide with red safety index mark 39 cm from distal tip



### Boyles Davis Mouth Gag (Set)

**Product description**

The Boyle Davis mouth gag with Tonsil plates, adult 5 blades. This instrument is part of the tonsillectomy set. Mouth gags are used to keep the patient's mouth open during oral surgery, leaving both hands of the surgeon free to operate. The Boyle Davis mouth gag consists of the Davis gag, a frame that serves to hold the mouth open and the Boyle tongue depressor to hold the tongue down. The tongue depressor comes in several sizes, from pediatric to adult. The instrument is assembled by sliding the tongue blade into the frame.



**Hey Groove Urethral Dilators Set of 4 Pcs.**



**Fine Tooth Forceps**



**Plain Fine Forceps**

**Optional Accessories**



**LED Hand held Light Source**

**Metzenbaum Scissor Curved 8 inch**





**Ryder Needle holder**  
4-0 8 inch



**Needle holder 3-0**  
6 inch T.C. Gold



**Hammer**



**Cheesle**



**Mosquito Forceps (5") with 1 - 6**  
marked for easy identification  
for suturing



**Perforated Box**

## Surgical Instrument for Radical Prostatectomy which Dr. S.K. Rawal is using in Rajiv Gandhi Cancer Institute, Delhi



**Apical Retractor**



**Prostate Retractor**



**Ryder Needle Holder**



**Prostatectomy  
Scissor Curved**



**Prostatectomy  
Scissor Half Curved**



**Needle Holder Curved (10")**



**Gemini Mixture Artery Forceps**



**Needle Holder Angular (10")**

**Manufacturer, Supplier, Importer & Exporter of :**

- UROLOGY ■ GENERAL SURGERY ■ GYNAECOLOGY ■ PLASTIC SURGERY
- DERMATOLOGY INSTRUMENTS ■ PAEDIATRIC SURGERY ■ E.N.T.
- OPHTHALMOLOGY ■ NEURO SURGERY ■ CARDIO THORACIS ■ LAPROSCOPY
- VETERINARY INSTRUMENTS & EQUIPMENTS ■ NURSERY EQUIPMENTS
- PHYSIOTHERPY EQUIPMENTS & SURGICAL COUTARY
- HOSPITAL FURNITURE & O.T. HOSPITAL EQUIPMENTS ETC.